



LG IMPLANTS TECHNICAL DATASHEET

Synthetic mesh for inguinal hernia by laparotomy



1. Administrative information about MICROVAL



MICROVAL
ZA CHAMP DE BERRE, 43240 SAINT JUST MALMONT, France
Tel: +33 4 77 35 03 03
Fax: +33 4 77 35 03 19

E-mail: info@microval.fr
Website: www.microval.fr



Medical device vigilance contact: Olivier CUILLERON
Tel: +33 4 77 35 03 03
Fax: +33 4 77 35 03 19
E-mail: info@microval.fr

2. Device information

2.1 Common name: Synthetic permanent implant for abdominal wall reinforcement

2.2 Commercial name: LG implant

2.3 Nomenclature code: GMDN 60300
EMDN P900202


2.4 Class of medical device: IIb according to European Directive 93/42/CEE (2007/47/CE)
Notified Body number: 1639
Date of first sale: 2007
Manufacturer: MICROVAL
Basic UDI-DI 37004584DT011-DIG-LAPWQ



Certificate No. : FR19-81843429

2.5 Device's description:








- ❖ These implants can be implanted by Lichtenstein procedure and its derivatives.
- ❖ Available in male or female versions (slot for spermatic cord or not).
- ❖ These permanent implants and their very high burst strength allow a durable inguinal reinforcement.

Characteristics	Value*
Mesh type	Knitted monofilament polypropylene – lock stitch
Construction	Standard
Thickness¹	0.56 mm
Weight²	90 g/m ²
Max pore size³	1.31 mm
Porosity⁴	≥77%
Burst resistance⁵ (max in vivo value ≈ 170mmHg ^a)	≥ 6135mmHg
Strain at Ultimate Tensile Test⁵	102% (longitudinal) 111% (transverse)
Durability	Permanent
Surgical technique	Laparotomy Lichtenstein
Microscopic view of the mesh	

*Average values given as an indication

¹ NF EN ISO 5084 (1996); ² ISO 3801 (1977) method 5
³ NF S94-801 (2007) method B; ⁴ NF S94-801(2007) method A
⁵ NF EN ISO 13938-1(2019); ⁵ NF EN ISO 13934-1 (2013)


^a Pott et al. 2012, « Mechanical Properties of Mesh Materials Used for Hernia Repair and Soft Tissue Augmentation », PLoS ONE 7(10): e46978. doi:10.1371/journal.pone.0046978

2.6	<div>References:</div> <table><tr><th>Size</th><th>Destination</th><th>Shape</th><th>Standard Polypropylene</th></tr><tr><td rowspan="2">6 x 12 cm</td><td>Woman ♀</td><td></td><td>413 135</td></tr><tr><td>Man ♂</td><td></td><td>413 136</td></tr></table>	Size	Destination	Shape	Standard Polypropylene	6 x 12 cm	Woman ♀		413 135	Man ♂		413 136
Size	Destination	Shape	Standard Polypropylene									
6 x 12 cm	Woman ♀		413 135									
	Man ♂		413 136									
2.7	<div>Device composition:</div> 100% polypropylene <div>✓ No latex</div> <div>✓ No phthalates</div> <div>✓ No products of animal or organic origin</div>											
2.8	<div>Field of use – Indications:</div> Abdominal wall reinforcement. Treatment of hernias and eventration.											
3. Sterilization												
3.1	<div>Sterilized:</div> <input checked="" type="checkbox"/> YES <input type="checkbox"/> NO											
3.2	<div>Sterilization process:</div> Ethylene Oxide according to NF EN ISO11135:2014 and NF EN ISO10993-7:2008(A1:2019)											
4. Storage conditions												
	<div>Packaging:</div> 1 implant packed in double Tyvek AND filmed cardboard box 240mm x 222mm x 20mm (non-contractual photograph) <div></div> <div>Expiration:</div> 5 years after sterilization <div>Storage:</div> no particular conditions, store at ambient temperature, please read D133 IFU											
5. Safety												
	Please read Instructions for Use D133											
6. Usage												
6.1	<div>IFU:</div> D133											
6.2	<div>Indication:</div> Abdominal wall reinforcement. Treatment of hernia and eventration.											
6.3	<div>Precautions of use:</div> Before operation, please check that all specific instruments for the operation are available and functional. Avoid any contact with objects which could damage the device. The damaged devices and/or that have been in contact with a patient must be isolated and disinfected before cleaning and possible back forwarding. Caution: a defect in the fixing or positioning of the device can induce abnormal stresses and/or reduce the service life.											
6.4	<div>Contra-indications:</div> Children during their growth, intensive and/or violent physical activities; Allergic reaction. Serious illness inducing a risk of dangerous post-operative complication. Infection and septicemia are absolute contra-indications.											
7. Additional information concerning the product												
	<div>Bibliography, test reports:</div> <div><div>❖ [Ref0510] “Standard polypropylene mesh vs lightweight mesh for Lichtenstein repair of primary inguinal hernia: a randomized controlled trial”, Z DEMTRASHVILI ,K KHUTSISHVILI, I PIPIA, G KENCHADZE, E EKALADZE, Int J of Surgery 2014, 12:1380-1384</div><div>❖ [Ref0509] “Open Mesh Techniques for inguinal hernia repair: a meta-analysis of randomized controlled trials”, G ZHAO, P GAO, B MA, J TIAN, K YANG, Annals of Surgery 2009, 250(1):35-42</div><div>❖ [Ref0020] “Etude rétrospective et analytique du traitement des hernies inguinales de l’adulte sur 130 patients de 1996 à 1997.”, J L DULUCQ and P WINTRINUGER, PMCF study Microval 1998</div></div>											

8. Appendices

8.1 IFU: D133

8.2 Labelling example:



8.3 Symbols used in IFU and/or labels:

