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(EN) REPOL ANGIMESH - ANGIMESH PRE - REPOL PLUG - FOLDED MESH - WINGS MESH

MESH FOR SURGERY INSTRUCTIONS FOR USE

DESCRIPTION, PERFORMANCE, GROUP OF PATIENTS / RECIPIENTS USERS

REPOL ANGIMESH, ANGIMESH PRE, REPOL PLUG (REPOL PLUG BASIC, REPOL PLUG CAP and REPOL PLUG CAP), FOLDED MESH and WINGS MESH surgical meshes are biocompatible prostheses made of 100% monofilament polypropylene intended for use in general surgery on adult patients suffering from abdominal hernial pathology (inguinal, femoral or crural and ventral), or for the repair / reinforcement / restoration of clinical conditions characterized by tissue laxity that require structural reinforcement through the implantation of a synthetic prosthetic material. All meshes are non-absorbable, resistant, easy to handle and, depending on the variant, they are available in various textures and sizes to better adapt to the needs of intervention.

Disposable, sterile, non-sterilizable devices. The devices do not contain latex, phthalates, products of animal or biological origin.

The devices must be used exclusively in the operating room, by surgeons experienced in the appropriate surgical techniques, familiar with the applicable international guidelines, as well as the relevant medical literature regarding surgical techniques, complications and dangers associated with the intended intervention procedures, in order to limit any unwanted effects and complications.

CLINICAL INDICATIONS AND CONDITIONS TO BE TREATED

REPOL ANGIMESH are meshes indicated for the treatment of abdominal hernias of various types, including the hernia, and for the reinforcement of soft tissues where there is the presence of a weak point. They are available in various shapes and sizes.

ANGIMESH PRE are pre-shaped meshes ideal for the surgical repair of inguinal hernia with anterior access. They are available both with a hole (to allow the passage of the spermatic cord) and without a hole.

REPOL PLUG, REPOL PLUG BASIC is ideal for the repair of inguinal hernias with anterior access. REPOL PLUG CAP and REPOL PLUG FLOWER are three-dimensional meshes with optimal softness, memory and consistency for the surgical repair of inguinal and femoral hernias with anterior access.

FOLDED MESH are pre-shaped meshes with self-locking flap that do not need to be closed with points of the final tails around the funiculus indicated for the treatment of inguinal hernias. This particular prosthesis can be used in association with a plug or separately depending on the technique that the surgeon wants to perform (sutures/ou tension free).

WINGS MESH are double layer meshes ideal for the surgical repair of inguinal hernias with anterior and posterior access. These meshes allow to perform the true "sutures/ou" technique as thanks to their wings it is possible to position them around the spermatic cord without the need for any fixation point.

The devices in the variants REPOL ANGIMESH, ANGIMESH PRE, REPOL PLUG BASIC and WINGS MESH can be implanted with the laparoscopic technique.

It is up to the surgeons to identify the most suitable technique to apply to the individual patient.

EXPECTED CLINICAL BENEFITS

The REPOL ANGIMESH, ANGIMESH PRE, REPOL PLUG (REPOL PLUG BASIC, REPOL PLUG CAP and REPOL PLUG CAP), FOLDED MESH and WINGS MESH surgical meshes, once implanted, allow the restoration of the integrity of the patient's damaged tissues, restoring the balance of normal performance of the same. By closing the hernial defect and forming a new tissue, these meshes allow the patient to resume normal life.

CONTRAINDICATIONS

Do not place the meshes in contact with the intestinal loops.

Do not use meshes on patients:

- with potential future growth or expansion of tissues (e.g., infants, children and women who are pregnant or considering future pregnancies);
- in potentially infected or clearly infected sites, or in the case of pre-existing acute or chronic infections;
- with pathologies of the soft tissues in which the device must be inserted;
- with any serious pathology that would compromise the patient's health and /or implant placement;
- subjected to anticoagulant therapy or with pathologies that limit blood supply and compromise scarring;
- with autoimmune diseases affecting connective tissues (for example Kohn's disease);
- with known allergy or sensitivity to the components.

SIDE EFFECTS

One of the side effects that could occur following the implantation of the device is the so-called foreign body reaction due to the consequent natural interaction between the material of which the device is made and the bio-chemical environment in which it is being implanted.

(FR) REPOL ANGIMESH - ANGIMESH PRE - REPOL PLUG - FOLDED MESH - WINGS MESH

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INFECTIONS that do not respond to antibiotic therapy may require partial or total removal of the prosthesis.

In some patients, the implantation of the device can cause the onset of postoperative pain. Other potential adverse reactions to device implantation are inflammation, bleeding and/or haemorrhage/hematoma, seroma formation, adhesion formation, fistula formation, extrusion or enhancement of infection, erosion and perforation of hollow viscera and organs, rejection of the device from the body.

To treat these adverse reactions, one or more revision surgeries may be required, and some complications may require removal of the device. Seromas are generally self-limiting and heal without the need for device removal.

The implantation of surgical meshes must always observe correct clinical practice of implantation and post-operative monitoring by the surgeon.

METHOD OF USE

Once the implant site has been identified, the mesh must be arranged in such a way as to adhere perfectly to the reference surfaces (muscular plane and fascial plane) without creating folds that reduce its contact surface. Thereafter, depending on the surgeon's technical preferences and the type of hernia, the mesh can be anchored with both resorbable and non-resorbable sutures.

In the case of implantation using the laparoscopic technique, trocars with dimensions between 5 and 12mm can be used.

WARNINGS AND PRECAUTIONS FOR USE (SEE ALSO CONTRAINDICATIONS)

Patients suitable for surgery and the choice of implant must be carefully selected after a complete preoperative diagnosis by the surgeon in order to limit any side effects and complications.

It is the surgeon's task to establish the best post-operative plan (treatments and check-ups) according to the specific case of each individual patient. In addition, the surgeon will have to take into consideration the patient's feedback on their quality of life and any undesirable effects encountered.

The device is supplied in a sterile condition and must be handled with care by opening its packaging in order to reduce the risk of contamination as much as possible.

During the surgery, care must be taken not to damage the large vessels, the bladder or other organs. The utmost attention must be paid to avoid that the site of implantation of the meshes is not exposed, or is near the skin surface, due to the risk of erosion or exposure and subsequent contamination.

ANGIOLOGICA cannot be held responsible in the event that one of the following cases occurs during the pre- and post-operative phase:

- insufficient experience or training of the surgeon;
- modification or adaptation of the device (for example a cutout);
- excessive traction which leads to a compromise of the physical and structural integrity of the device;
- surgical technique not performed correctly;
- lack of knowledge of the applicable international guidelines, of the commercial and / or technical documents of the device and instruments;
- unsuitable therapies or insufficient explanation of post-operative precautions;
- failure to comply with the contents of these instructions.

Before use, the device is comparable to solid urban waste. If an unused device comes into contact with instruments or other surgical materials used on a patient, or if it has been contaminated with organic fluids, discard it carefully to prevent the risk of transmitting a viral infection. After use, the medical device must be disposed of as sanitary waste and in case of explant dispose of the system in accordance with the procedure applicable to biological waste. The device can be stored at room temperature and the general storage conditions must be suitable for maintaining the integrity of the device itself. This device complies with the manufacturer's requirements if delivered and properly stored in the original packaging. ANGIOLOGICA assumes no responsibility in the event that changes are made to the original packaging.

Do not use after the expiry date indicated on the package.

The device does not affect and is not affected by magnetic resonance imaging (MRI).

United States federal law permits this product to be sold only to or on the order of physicians.

It is necessary to report any serious accident occurring in relation to the medical device supplied to the manufacturer and to the competent authority of the Member State in which it is based.

STERILIZATION

Ethylene Oxide (EO) sterilized device.

The contents are sterile if the package has not been opened or damaged. Do not use the device if the package is opened or damaged or has been accidentally opened before use.

The product is for single use only and must not be re-sterilized as:

- for use on a single patient;
- not reusable or reesterilizable, since otherwise:
 - the structural integrity of the device may be compromised and it could cause failure, with a consequent risk of injury;
 - there is a risk of device contamination and/or patient infections and cross-infection, including the transmission of infectious diseases.

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SSCP (SUMMARY OF SAFETY AND ANCLINAL PERFORMANCE)

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