(IT) COMBI MESH PLUS - COMBI MESH PLUS PRO

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MESH FOR SURGERY INSTRUCTIONS FOR USE

DESCRIPTION, PERFORMANCE, GROUP OF PATIENTS / RECIPIENTS USERS

COMBI MESH PLUS and COMBI MESH PLUS PRO surgical meshes are biocompatible double-surface prostheses made of 100% monofilament polypropylene (PP) coated with polyurethane (PU) intended for use in general surgery on adult patients suffering from abdominal hernial pathology (inguinal, femoral or crural and ventral), or for the repair / reinforcement / restoration of clinical conditions characterized by tissue laxity that require structural reinforcement through the implantation of a synthetic prosthetic material. All meshes are non-absorbable, resistant, easy to handle and are available in various sizes to better adapt to the needs of intervention.

Disposable, sterile, non-sterilizable devices. The devices do not contain latex, phthalates, products of animal or biological origin. The devices should only be used in the operating room, by surgeons experienced in the appropriate surgical techniques, familiar with the applicable international guidelines, as well as the reference medical literature regarding surgical techniques, complications and dangers associated with the planned intervention procedures.

CLINICAL INDICATIONS AND CONDITIONS TO BE TREATED

COMBI MESH PLUS is indicated for the repair of major defects of the abdominal wall especially in the treatment of laparocels or in the closure of the abdominal wall in all those cases in which the integrity of the wall is compromised.

COMBI MESH PLUS PRO are the ideal meshes for surgical repair of inguinal hernia with anterior access. They are already ready for use as their particular shape with perfectly finished edges allows the surgeon to avoid an intraperitic shaping of the prosthesis. The hole allows the passage of the spermatic cord without creating adhesions. It is up to the surgeons to identify the most suitable technique to apply to the individual patient.

EXPECTED CLINICAL BENEFITS

COMBI MESH PLUS and COMBI MESH PLUS PRO surgical meshes, once implanted, allow the restoration of the integrity of the patient's damaged tissues, restoring the balance of normal performance of the same. By closing the hernial defect and forming a new tissue, these networks allow the patient to resume normal life.

CONTRAINDICATIONS

Do not place the PP surface in contact with the intestinal loops.

Do not use meshes on patients:

- with potential future growth or expansion of tissues (e.g., infants, children and women who are pregnant or considering future pregnancies);
- in potentially infected or clearly infected sites, or in the case of pre-existing acute or chronic infections;
- with pathologies of the soft tissues in which the device must be inserted;
- with any serious pathology that would compromise the patient's health and /or implant placement;
- subjected to anticoagulant therapy or with pathologies that limit blood supply and compromise scarring;
- with autoimmune diseases affecting connective tissues (for example Kohn's disease);
- with known allergy or sensitivity to the components.

SIDE EFFECTS

One of the side effects that could occur following the implantation of the device is the so-called foreign body reaction due to the consequent natural interaction between the material of which the device is made and the bio-chemical environment in which it is being implanted.

Infections that do not respond to the antibiotic therapy may require the partial or total removal of the prosthesis.

In some patients, the implantation of the device can cause the onset of postoperative pain.

Other potential adverse reactions to device implantation are inflammation, bleeding and/or haemorrhage/hematoma, seroma formation, adhesion formation, fistula formation, extension or enhancement of infection, erosion and perforation of hollow viscera and organs, rejection of the device from the body. To treat these adverse reactions, one or more revision surgeries may be required, and some complications may require removal of the device. Seromas are generally self-limiting and heal without the need for device removal.

The implantation of surgical meshes must always observe correct clinical practice of implantation and post-operative monitoring by the surgeon.

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METHOD OF USE

Once the implant site has been identified, the mesh must be arranged in such a way as to adhere perfectly to the reference surfaces (muscular plane and fascial plane) without creating folds that reduce its contact surface. During surgery, the COMBI MESH PLUS mesh should be positioned with the PP surface facing the wound opening and the PU surface towards the inside of the patient's body. Instead, the COMBI MESH PLUS PRO mesh should be positioned with the PP surface facing inwards of the patient's body and the PU surface towards the wound opening and in contact with the spermatic cord.

Later, depending on the clinical preferences of the surgeon and the type of hernia, the mesh can be anchored with both resorbable and non-absorbable sutures. In case of implantation by laparoscopic technique, trocars with dimensions between 5 and 12 mm can be used.

WARNINGS AND PRECAUTIONS FOR USE (SEE ALSO CONTRAINDICATIONS)

To help the correct positioning of the mesh, a colored easily removable PP wire with visible ends has been introduced on the PP side.

Patients suitable for surgery and the choice of implant must be carefully selected after a complete preoperative diagnosis by the surgeon in order to limit any side effects and complications.

It is the surgeon's task to establish the best post-operative plan (treatments and check-ups) according to the specific case of each individual patient. In addition, the surgeon will have to take into consideration the patient's feedback on their quality of life and any undesirable effects encountered.

The device is supplied in a sterile condition and must be handled with care by opening its packaging in order to reduce the risk of contamination as much as possible.

During the surgery, care must be taken not to damage the large vessels, the bladder or other organs.

The utmost attention must be paid to avoid that the site of implantation of the network is not exposed, or is near the skin surface, due to the risk of erosion or exposure and subsequent contamination.

ANGIOLOGICA cannot be held responsible in the event that one of the following cases occurs during the pre- and post-operative phase:

- insufficient experience or training of the surgeon;
- modification or adaptation of the device (for example a cutout);
- excessive traction which leads to a compromise of the physical and structural integrity of the device;
- surgical technique not performed correctly;
- lack of knowledge of the applicable international guidelines, of the commercial and / or technical documents of the device and instruments;
- unsuitable therapies or insufficient explanation of post-operative precautions;
- failure to comply with the contents of these instructions.

Before use, the device is comparable to solid urban waste. If an unused device comes into contact with instruments or other surgical materials used on a patient, or if it has been contaminated with organic fluids, discard it carefully to prevent the risk of transmitting a viral infection. After use, the medical device must be disposed of as sanitary waste and in case of explant disposal of the system in accordance with the procedure applicable to biological waste.

The device can be stored at room temperature and the general storage conditions must be suitable for maintaining the integrity of the device itself. This device complies with the manufacturer's requirements if delivered and properly stored in the original packaging. ANGIOLOGICA assumes no responsibility in the event that changes are made to the original packaging.

Do not use after the expiry date indicated on the package.

The device does not affect and is not affected by magnetic resonance imaging (MRI).

United States federal law permits this product to be sold only to or on the order of physicians.

It is necessary to report any serious accident occurring in relation to the medical device supplied to the manufacturer and to the competent authority of the Member State in which it is based.

STERILIZATION

Ethylene Oxide (EO) sterilized device. The contents are sterile if the package has not been opened or damaged. Do not use the device if the package is opened or damaged or has been accidentally opened before use.

The product is for single use only and must not be re-sterilized:

- for use on a single patient;

- not reusable or re-sterilizable, since otherwise:
 - the structural integrity of the device may be compromised and/or its failure could be caused, with a consequent risk of injury;
 - there is a risk of device contamination and/or patient infections and cross-infection, including the transmission of infectious diseases.

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