Fabricant : MICROVAL
ZA Champ de Berre – 43240 Saint-Just-Malmont – France
Tel : +33(0)4 77 35 03 03 Fax : +33(0)4 77 35 03 09
e-mail: info@microval.fr www.microval.fr

FABRIQUÉ
EN FRANCE

MICROVAL France
Conception et fabrication de dispositifs médicaux –
Medical devices design and manufacturing

MD

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
* see recharging methods in paragraph 12.

DESCRIPTION

SPIRE'IT® is a device comprising of an applicator and sutures for fitting wall reinforcement implants as part of treatment for hernias or prolapses (abdominal approach). The SPIRE'IT® applicator allows the implantation up to 20 sutures in laparoscopic and laparatomic surgery. The sutures are made of a nickel titanium shape memory alloy (NITINOL) in the form of double 4mm in diameter coils. The MICROVAL SPIRE'IT® applicator is a sterile device designed for multiple use during a single procedure. The long SPIRE'IT® applicators for laparoscopic use are to be used using a trocar with a minimum diameter of 5mm. The MICROVAL SPIRE'IT® applicator is delivered ready to use, with 10 fixing coils and an additional 10 coil charger (for the references concerned). The applicator implants a suture each time the trigger is pressed (A). On leaving its casing, the suture closes naturally to form a coil, fixing the implant in the tissues. (A) trigger, (B) end, (C) magazine inside B, (D) rotary indicator on back of applicator. The applicator is a device for temporary use; sutures are permanent implants.

1. USERS WHO THE DEVICE IS INTENDED FOR

The use of MICROVAL devices is reserved for practitioners specializing* in the fitting of reinforcement implants to treat inguinal hernias, in order to limit undesirable effects and complications.

This device must only be used in suitable, aseptic premises.

* Practitioners who have completed sufficient training to fit these implants via PET or TAPP and who practice these procedures regularly

2. TARGET PATIENT GROUPS

The device is intended for adults with weakness of the abdominal wall such as a hernia or eventration and for adult women with advanced apical pelvic organ prolapse (POP-Q ≥ stage 2) in the setting of abdominal route surgery.

3. INDICATIONS AND MEDICAL CONDITIONS TO BE TREATED

The device is indicated in the context of implant fitting for laparoscopic or laparatomic hernia repair.

The device is indicated in the context of implant fitting for treatment of abdominal prolapse by laparoscopic or laparatomic route.

4. CONTRAINDICATIONS

SPIRE'IT® cannot be used alone; its contraindications are those regarding fitting of a wall reinforcement implant:

- wall reinforcement implants for hernia and eventration repair are not indicated for pregnant women, in growing children, in adults with intense physical activity, in adults who are allergic to nitinol, in adults with serious diseases with a risk of post-operative complications, and in adults with infection or sepsis.
- wall reinforcement implants to treat prolapse (abdominal route) are not indicated for pregnant women or women planning a future pregnancy, in growing children, in women with intense physical activity, in women who are allergic to nitinol, in women with serious illnesses with a risk of post-operative complications, and in women with infection or sepsis.

Do not use the SPIRE'IT® applicator on vascular or nervous tissues, or on solid organs (bones). However, it is possible to place fixations above a bone, underlying vessels or viscera, but in this case, there should be a gap greater than 4.5 mm between the tissue surface and the underlying structure.

Do not use the SPIRE'IT® applicator to achieve haemostasis.

Do not use the SPIRE'IT® applicator on a site where it is impossible to visually observe haemostasis.

Avoid prolonged exposure to excessive temperatures.

5. PRECAUTIONS

5.1. Precautions for use

The surgeon must check, before the operation, that the specific instruments for the use of the devices are available and functional.

SPIRE'IT® is delivered sterile (sterilized by GAMMA radiation). Avoid any contact with objects that could damage the device (applicator and sutures). If the device or its sterile packaging is damaged before the procedure, the device must not be fitted.

Warning: a lack of fixation, positioning of the implant, a lack of use of SPIRE'IT® can induce abnormal stress and/or reduce the useful life of the device and/or reduce the effectiveness of the wall reinforcement implant.

The peri- and post-operative results obtained by the user and the medical team may not be appropriate in cases in which there is:

- insufficient practitioner experience and/or training ;
- a modification or transformation of the device (cutting for example) ;
- excessive traction resulting in degradation of the integrity of the device ;
- incorrect indication ;
- faulty surgical technique ;
- a lack of knowledge of the standards in force, of the commercial and/or technical documents, of the system, of the instruments ;
- a lack of care or a lack of rehabilitation after the operation ;
- failure to comply with the contents of this manual.

In these cases, the user of this device will be held solely responsible, thus relieving MICROVAL of any liability.

Applying excessive pressure to the end of the pliers while in use can cause it to jam. After any suture coil implantation using the pliers, always carefully check haemostasis. It is then possible to treat bleeding by electrocoagulation or by manual suture. It is not advisable to reload the SPIRE'IT® applicator more than once during the same procedure, otherwise the integrity of the fixation will not be guaranteed. The SPIRE'IT® applicator is designed for multiple use during a single procedure.

5.2. Side effects and complications

Seroma/haematoma formation; urinary retention; neuralgia; testicular pain; oedema; infection (implant or surgical site); recurrent hernia, and eventration or prolapse.

5.3. Residual risks

The residual risks are the risks that persist despite all the measures implemented by MICROVAL, these are the risks associated with the user of the device:

- The medical team modifies the device ;
- The implant is incorrectly positioned and causes pain to the patient ;
- Instruments and consumables used in the operating room are no longer sterile ;
- The implant is contaminated during the operation ;
- The surgeon misjudges the general level of health of the patient ;
- The surgeon's surgical procedure and/or diagnosis lack precision.

Risks related to the environment of the device:

- The operating environment is not sterile ;
- The temperature during storage or during transport is too low or too high and may influence the conformity of the product.

Risks related to non-conformities not detected at MICROVAL.

6. CLINICAL BENEFIT AND PERFORMANCE

The performances and the expected clinical benefits for the patient, claimed by MICROVAL are an optimal cure rate and patient satisfaction rate, a rapid return to normal activity, and the shortest possible hospital stay.

The performances and the clinical benefits expected for the user, claimed by MICROVAL are the shortest possible intervention duration and the lowest possible pain in the short and medium term.

7. INSTRUCTIONS NEEDED TO HELP THE USER TO USE THE DEVICE CORRECTLY

7.1. Preoperative evaluation of the patient

The intervention must be preceded by an initial evaluation in consultation. The decision to perform a medical act must be taken in agreement with the duly informed patient, who has had sufficient time to reflect. The user must carefully consider which implant to use, taking into account the indication encountered, the morphology, the activity, the physical and/or mental state of the patient, to limit unwanted side effects and complications. SPIRE'IT® (20) offers the possibility of placing 20 sutures for the fixation of an implant; it is the surgeon's responsibility to ensure that the stitches are sufficient, while at the same time using the least amount of implanted material.

7.2. Surgical intervention

An explanatory video is available at www.microval.fr.

Check that all the instruments are compatible with each other. Remove the SPIRE'IT® applicator from its packaging using the usual techniques for sterile products. Do not tilt the applicator in the sterile field to avoid damaging it. In the case of laparoscopic use, insert the applicator into the appropriate trocar (minimum diameter 5mm). Position the tip of the applicator centred over the surgical site and perpendicular to the implant.

Apply light pressure to the tissue site while the suture is being placed. When placing an internal suture, if possible, apply external back pressure to the area to be fixed to ensure a more secure fixation.

To actuate the applicator, press the tip against the tissue and/or the implant and at the same time squeeze the trigger (A) fully, then release, which prepares the applicator for the next application. Do not suddenly release the trigger while applying a suture. Do not apply excessive force when triggering the applicator as it may jam. If a suture jams, remove the entire instrument. Examine the stapled area to ensure the correct formation and position of the coil. Repeat all the above steps over the entire surface of the implant to be stapled until the operation is completed.

SINGLE USE CHARGER (C):

A rotary indicator (D) located at the back of the applicator allows you to locate the status of the charger. At the start of the use of the applicator, the initial position of this indicator must be noted. When this indicator has made a complete revolution, the 10 fixing turns have been installed, it is then possible to change the magazine.

The two laparoscopic or laparotomic methods recommended for the placement of wall reinforcement implants are used at the discretion of the surgeon according to the recommendations of the instruction leaflets relating to these implants. For MICROVAL implants, refer to instructions D121, D133 and D167 for digestive implants. MICROVAL provides a description of the use of SPIRE'IT® on request.

7.3. Post-operative evaluation

A follow-up consultation should be performed within one month of implantation. During this consultation, patients' feedback on their quality of life and the adverse events experienced must be considered in order to detect and deal with any complications early on. MICROVAL recommends a consultation one year after the operation, in order to ensure active management of any late complications. Other consultations can be carried out as needed.

The management of serious post-implantation complications is the subject of multidisciplinary consultation and a decision shared with the patient after they have been duly informed of all the options and has benefited from a period of sufficient reflection. If an explant is necessary, it must be performed in a centre with a multidisciplinary surgical technical team and must be reserved for surgeons trained in the explantation of implants. This training is based on practical training acquired by learning from an experienced surgeon.

8. MATERIALS AND SUBSTANCES TO WHICH PATIENTS ARE EXPOSED

Absence of latex, phthalates, animal or biological products.

9. STORAGE / HANDLING / DISPOSAL

This device must be delivered and stored in the original MICROVAL packaging and stored at room temperature. Storage conditions must avoid all risks of alteration of the labelling, packaging, and device. The packaging must not be modified or transformed. In the event of an explant, the implant should be disposed of in accordance with the procedure for biological waste. The accessories with which the device is used, if they cannot be re-sterilized, must be disposed of in the same way.

For safety reasons, any applicator to be returned to MICROVAL must be decontaminated if it has been soiled.

10. INFORMATION TO BE GIVEN TO THE PATIENT IF THEY ARE NOT THE USER

The user of this device must be aware of and must inform the patient before the procedure:


- of the precautions to be taken after the operation as well as the contraindications and undesirable events that may occur ;
- that this device may only allow a lower level of activity, certain activities may be partially or totally prohibited over time (defined by the user) ;
- that the inevitable wear and tear of the device, possible complications, and/or possible detachment over time may force the user or another surgeon to intervene again.

In the event of material vigilance, or when a serious incident occurs, the surgeon or the patient must notify the manufacturer and/or the authorities as soon as possible.

11. INFORMATION RELATING TO THE REUSE OF THE MEDICAL DEVICE

This device is supplied sterile (sterilized by GAMMA radiation) and for single use. It should not be reused or re-sterilized. Once the sterile packaging is opened, MICROVAL can no longer guarantee the sterility of the product; the implant then poses a risk of contamination for the patient.

12. CARTRIDGE CHANGE PROCEDURES (only for references with additional 10 cartridges, i.e., 943052 and 943302)

CAUTION:

⇒ CHECK THE ALIGNMENT OF THE ARROWS. The applicator is ready to operate again.

⇒ It is not appropriate to recharge the applicator more than once for the same procedure. Otherwise, the clamp could no longer function properly and jam or lead to poor fixation.

LABEL SYMBOLS:

STERILE R	GAMMA radiation sterilization		Do not reuse		Manufacturer
LOT	Batch number		Do not re-sterilize		Manufacture date YYYY-MM
REF	Product reference		Do not use if the packaging is damaged		Applicator length
	Caution		Expiry date		Number of coils available
	Refer to Instruction Manual	MD	Medical device		
UDI	Unique Device Identifier		Unique sterile barrier system		

	Hersteller: MICROVAL ZA Champ de Berre - 43240 Saint-Just-Malmont – Frankreich Tel.: +33 (0) 4 77 35 03 03 Fax: +33 (0) 4 77 35 03 09 E-Mail: info@microval.fr www.microval.fr	FABRIQUÉ EN FRANCE 	MICROVAL France <small>Conception et fabrication de dispositifs médicaux – Medical devices design and manufacturing</small>
MD	Gerät: SPIRE'IT® Nahtapplikator		D132v7 (DE)

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