	Fabricant : MICROVAL ZA Champ de Berre – 43240 Saint-Just-Malmont – France Tel : +33(0)4 77 35 03 03 Fax : +33(0)4 77 35 03 09 e-mail: info@microval.fr www.microval.fr	FABRIQUÉ EN FRANCE 	
	Dispositifs MICROVAL 2D MESH / 2F MESH	D167v0.1	(FR)

DESCRIPTION

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DESCRIPTION

MICROVAL 2D MESH / 2F MESH devices are wall reinforcement implants that mechanically reinforce the patient's tissues. These implants are made of polypropylene or knitted polyester. They can be rolled up to be introduced into a trocar. These implants retain a shape memory, allowing easy unrolling and perfect adaptation to the abdominal wall. To facilitate implantation, a black biocompatible silicone ink mark can be affixed to the implant. The devices are permanent implants.

1. USERS WHO THE DEVICE IS INTENDED FOR

The use of MICROVAL devices is reserved for practitioners specializing* in the fitting of reinforcement implants to treat all types of hernias (inguinal hernias, umbilical hernias, femoral hernias, eventrations) in order to limit undesirable effects and complications. This device must only be used in suitable, aseptic premises.

** surgeons experienced in the appropriate surgical techniques, familiar with the applicable international guidelines, as well as the relevant medical literature regarding surgical techniques, complications and dangers associated with the intended intervention procedures, in order to limit any unwanted effects and complications.*

2. TARGET PATIENT GROUPS

The device is intended for adults with weakness of the abdominal wall such as a hernia or an eventration.

3. INDICATION AND MEDICAL CONDITIONS TO TREAT

The device is indicated in the treatment of any type of hernia: inguinal hernia, crural or femoral hernia, umbilical hernia, eventration, and is fitted via laparoscopic (TAPP or TEP) or laparotomic route.

4. CONTRAINDICATIONS

The device is contraindicated:

- In pregnant women,
- In growing children,
- In patients with intense physical activity,
- In patients who are allergic to polypropylene and/or polyester,
- In patients with serious diseases posing a risk of dangerous postoperative complications,
- In patients with infection or sepsis.

Previous surgery on the lower abdomen is often considered a relative contraindication to the PET technique, as access to the preperitoneal space may be difficult if scar tissue is present from previous operations. The surgeon may then consider repairing an open inguinal hernia with another type of implant.

5. PRECAUTIONS

5.1. Precautions for use

The use of the TEP technique over the TAPP method should be limited to specialists/surgeons with appropriate training. Surgeons must be able to advise and inform their patients in an appropriate manner about the different treatments and the benefit/risk ratio of using the implant, before making the decision to undergo surgery.

Patients with inguinal hernia should undergo laparoscopic repair, unless there is reason to believe that an open repair might be more clinically appropriate (see previous paragraph regarding contraindications for example).

Surgeons must check, before the operation, that the specific instruments for fitting the device are available and functional.

The device is delivered sterile; the sterilization method (ethylene oxide for polypropylene implants or GAMMA radiation for polyester implants) is indicated on the label.

Avoid any contact with objects that could damage the device.

If the device or its sterile packaging is damaged before the procedure, the device must not be fitted.

The peri- and post-operative results obtained by the user and the medical team may not be appropriate in cases in which there is:

- insufficient practitioner experience and/or training,
- a modification or transformation of the device (cutting for example),
- excessive traction resulting in degradation of the integrity of the device,
- incorrect indication,
- faulty surgical technique,
- a lack of knowledge of the standards in force, of the commercial and/or technical documents, of the system, or of the instruments,
- a lack of care or a lack of rehabilitation after the operation,
- failure to comply with the contents of this manual.

In these cases, the user of this device will be held solely responsible, thus relieving MICROVAL of any liability.

Warning: a defect in fixing and/or fitting the device can cause abnormal stress and/or reduce its lifespan and/or reduce its effectiveness.

5.2. Adverse side effects and complications

Hematoma formation, delayed healing, sensitization reaction to implanted materials, device change in behavior or device migration, surgical site infection, mesh infection, adhesions, seroma formation, paresthesia, neuralgia, urinary retention, intestinal injury, intestinal/bowel obstruction, fever, bleeding, hydrocele, nerve damage, orchitis, acute and/or chronic pain, male infertility, testicular atrophy, endourethral resection of the prostate, appearance of another hernia, recurrence of hernia, hemorrhage, discomfort, Injury (bladder, intestinal, vas deferens, spermatic cord, vessels, nerve), hypoesthesia, testicular pain, emotional distress, death .

5.3. Residual risks

The residual risks are the risks that persist despite all the measures implemented by MICROVAL, these are the risks associated with the user of the device:

- The medical team modifies the implant,
- The implant is incorrectly fitted and causes pain to the patient,
- Instruments and consumables used in the operating room are no longer sterile,
- The implant is contaminated during the operation,
- The surgeon misjudges the general level of health of the patient,
- The surgeon's surgical procedure and/or diagnosis lack precision.

Risks related to the environment of the device:

- The operating environment is not sterile,
- The temperature during storage or during transport is too low or too high and may influence the conformity of the product.

Risks related to non-conformities not detected at MICROVAL.

6. CLINICAL BENEFIT AND PERFORMANCE

The performances and the clinical benefits expected and claimed by MICROVAL are a low rate of recurrence, a low rate of chronic pain and an optimal rate of patient satisfaction.

7. INSTRUCTIONS NEEDED TO HELP THE USER TO USE THE DEVICE CORRECTLY

7.1. Preoperative patient assessment

The intervention must be preceded by an initial evaluation in consultation.

The decision to fit a wall reinforcement implant must be taken in agreement with the duly informed patient, who has had sufficient time to reflect.

The user must carefully consider which implant to use, taking into account the indication encountered, the morphology, the activity, the physical and/or mental state of the patient, to limit unwanted side effects and complications.

7.2. Surgical intervention

Surgical techniques are used at the discretion of the surgeon in relation to the MICROVAL recommendations described in Table 1 above, for the fitting of implants. If MICROVAL SPIRE'IT instruments are used for the fitting of these implants, please refer to the SPIRE'IT instruction manual (D132).

7.3. Post-operative evaluation

At the end of the intervention, the implant card is given to the patient for their information and medical device traceability.

A follow-up consultation should be performed within one month of implantation. During this consultation, patients' feedback on their quality of life and the adverse events experienced must be considered in order to detect and deal with any complications early on. MICROVAL recommends a consultation one year after the operation, in order to ensure active management of any late complications. Other consultations can be carried out as needed.

The management of serious post-implantation complications is the subject of multidisciplinary consultation and a decision shared with the patient after he has been duly informed of all the options and has benefited from a period of sufficient reflection. If an explant is necessary, it must be performed in a centre with a multidisciplinary surgical technical team and must be reserved for surgeons trained in the explantation of wall reinforcement implants. This training is based on practical training acquired by learning from an experienced surgeon.

8. MATERIALS AND SUBSTANCES TO WHICH PATIENTS ARE EXPOSED

100% polypropylene or 100% polyester / Black biocompatible silicone ink / No latex, phthalates, animal or biological products.

9. STORAGE / HANDLING / DISPOSAL

This device must be delivered and stored in the MICROVAL packaging, the storage conditions (ambient temperature) must avoid all risks of altering the labelling, the packaging and the device. The packaging must not be modified or transformed.

Damaged devices and/or those that have been in contact with a patient must be isolated and decontaminated before cleaning and possible return to MICROVAL.

In the event of an explant, the implant should be disposed of in accordance with the procedure for biological waste. The accessories with which the device is used, if they cannot be re-sterilized, must be disposed of in the same way.

10. INFORMATION TO BE GIVEN TO THE PATIENT IF THEY ARE NOT THE USER

The user of this device must be aware of and must inform the patient before the procedure:



- of the precautions to be taken after the operation as well as the contraindications and undesirable events that may occur,
- that this device may only allow a lower level of activity, certain activities may be partially or totally prohibited over time (defined by the user),
- that the inevitable wear and tear of the device, possible complications, and/or possible detachment over time may force the user or another surgeon to intervene again.

In the event of material vigilance, or when a serious incident occurs, the surgeon or the patient must notify the manufacturer and/or the authorities as soon as possible.

11. INFORMATION RELATING TO THE REUSE OF THE MEDICAL DEVICE

This device is supplied sterile (sterilized with ethylene oxide for polypropylene implants and GAMMA radiation for polyester implants) and for single use. It should not be reused or re-sterilized. Once the sterile packaging is opened, MICROVAL can no longer guarantee the sterility of the product; the implant then poses a risk of contamination for the patient.

LABEL SYMBOLS:

	Sterilization by means of ethylene oxide		Do not re-sterilize		Standard knitted polypropylene (PPT) / lightweight mesh (PPT-LW)
	GAMMA radiation sterilization		Do not use if the packaging is damaged		Knitted polyester (PET) / Three-dimensional (PET-TRIDI)
	Batch number		Expiry date		Thermoformed standard polypropylene (BULMESH / BUMPMESH)
	Product reference		Medical device		Unique Device Identifier
	Caution		Double sterile barrier system		
	Refer to Instruction Manual		Manufacturer		
	Do not reuse		Manufacture date		
	Store out of direct sunlight		Temperature limit (°C)		
	Resists humidity				

	Hersteller: MICROVAL - ZA Champ de Berre - 43240 Saint-Just-Malmont - Frankreich Tel.: +33(0)477350303 / Fax: +33(0)477350309 E-Mail: info@microval.fr www.microval.fr	FABRIQUÉ EN FRANCE 	
	Produkte MICROVAL 2D MESH / 2F MESH	D167v0.1	(DE)

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