MICROVAL

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ENGLISH / US

DESCRIPTION

The Disposable Suture Needle is composed of a plunger and needle with suture grasper.

INDICATION FOR USE

The Suture Grasper Needle is used to facilitate the placement of sutures during the closure of laparoscopic wound sites and laparoscopic fixation of mesh materials used in the repair of hernia defects or soft tissue deficiencies.

CONTRAINDICATIONS

The device is not intended for use when laparoscopic and other minimally invasive surgical procedures are contraindicated.

CAUTION

For single patient use only. The design of this device may not perform as intended by the manufacturer if it is re-used. The manufacturer cannot guarantee the performance, safety and reliability of a reprocessed device. After use, this product may be a potential biohazard. Handle in a manner, which will prevent accidental puncture. Dispose of in accordance with applicable laws and regulations. Carefully place the used needle in a sharps biohazard container after the procedure is completed.

WARNINGS AND PRECAUTIONS

The Disposable Suture Grasper Needle is provided in a sterile packaging. The packaging is to be checked for any damage before use. Do not use the device if the packaging is open or damaged. Extreme care should be used during insertion to avoid inadvertent puncture of any internal organs. Procedure should be performed only by qualified and trained physicians familiar with laparoscopic surgical techniques.

DIRECTIONS OF USE

Placement of sutures during the closure of laparoscopic wound sites

1. Acquire a suture of sufficient length. Remove the needle guard from the Suture Grasper Needle. The arrow indicator mark is aligned with and designates the beveled side of the needle. Depress the plunger on the proximal end of the device to advance the grasper arms past the distal end of the needle. Place the midpoint of the suture between the grasping arms. Release the plunger to retract the suture into the needle. Pull the two lengths of suture back along the needle towards the handle end of the device and hold in place with gentle tension.
2. With the trocar still in place, position the needle at one edge of the wound site. Angle the trocar in a position that allows for precise needle advancement. Hold the needle at the preferred angle and advance the needle through the subcutaneous tissue and fascia. Keep the needlepoint in direct visualization after it has entered the abdominal cavity. Briefly depress the plunger to release the suture. Withdraw the device.
3. Position the needle at the opposite edge of the trocar wound. Angle the trocar in a position that allows for precise needle advancement. Hold the needle at the preferred angle and advance the needle through the subcutaneous tissue and fascia. Keep the needlepoint in direct visualization after it has entered the abdominal cavity. Depress the plunger to advance the grasping arms. Grasp the previously placed suture loop with the arms. Release the plunger to retract the suture into the needle. Hold one of the two ends of the suture to prevent movement. Withdraw the needle to pull the suture through the tissue. Continue to retract the suture loop until the other end of the suture exits the grasping arms.
4. Remove the trocar from the wound site. Using manual traction on the suture withdraw any excess length that may remain in the abdominal cavity. Use standard technique to tie the suture.
5. Dispose in accordance with applicable laws and regulation. Carefully place the used needle in a sharps biohazard container after the procedure has been completed.

Transparietal fixation of mesh materials

1. Before the introduction of the mesh, position sutures on the wall side of the mesh material. The length of both arms of the suture should be adapted to the thickness of the abdominal wall.
2. Remove the needle guard from the Suture Grasper Needle. The arrow indicator mark is aligned with and designates the beveled side of the needle.
3. When the mesh is well positioned inside the abdominal cavity, advance the Suture Grasper needle through the abdominal tissue and into the body cavity under laparoscopic visualization.
4. Depress the plunger on the proximal end of the device to advance the grasper arms past the distal end of the needle.
5. Place the midpoint of the suture between the grasping arms. Release the plunger to retract the suture into the needle.
6. Withdraw the Suture Grasper Needle.
7. Depress the plunger to release the suture. Grasp and retain the suture with an instrument.
8. Repeat steps 2-6 for placement of the second arm of the suture as desired.
9. Tie the suture.
10. Dispose in accordance with applicable laws and regulation. Carefully place the used needle in a sharps biohazard container after the procedure has been completed.

PACKAGING

Packed individually in a sterile pouch. Provided in a box of 10 units. Single use device. Sterilized by ethylene oxide.

STORAGE

Do not use the device past the last day of the labeled month of expiration. Upon receipt of device, ensure that the packaging is not open or damaged and retains its sealed integrity. Do not use the device if the package appears to be damaged in any way that could affect sterility.

CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

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CE 0086  Ranfac Corp.