



# INGUINAL HERNIA CŒLIOSCOPY



Precut implant for TEP  
or TAPP surgeries.

90G / M<sup>2</sup>  
POROSITY > 77%



60G / M<sup>2</sup>  
POROSITY > 84%



100G / M<sup>2</sup>  
POROSITY > 80%



	cm	STANDARD POLYPROPYLENE	LIGHTWEIGHT POLYPROPYLENE	STANDARD POLYESTER
	10x14	413 410	414 410	416 410
	11x15	413 511	414 511	416 511
	12x16	413 612	414 612	416 612