




INGUINAL HERNIA

CÆLIOSCOPY



Tridimensional implant
for TEP or TAPP surgeries.

90G / M²
POROSITY > 77%



60G / M²
POROSITY > 84%



PICTURES	SIZES	STANDARD POLYPROPYLENE	LIGHTWEIGHT POLYPROPYLENE
	10x14	411 040 ^R 411 041 ^L 411 140 ^S	412 040 ^R 412 041 ^L 412 140 ^S
	11x15	411 050 ^R 411 051 ^L 411 150 ^S	412 050 ^R 412 051 ^L 412 150 ^S
	12x16	411 060 ^R 411 061 ^L 411 160 ^S	412 060 ^R 412 061 ^L 412 160 ^S

R - right / L - left / S - symmetrical